

S.P.I.N.E Plus

Code L0627 Approved

Features and Benefits:

- Patented, new posterior panel and SPINE brace that comfortably contour and conform to each patient's body structure.
- 6:1 ratio, low friction Mechanical Advantage Pulley System provides powerful, smooth and easily-controlled compression.
- Lightweight, breathable construction combined with contouring panels for increased comfort and support.
- Quick, one-hand adjustment controls support level and comfort for standing or sitting.
- Single front closure provides easy opening and closing.
- 2 panel sizes available.
- New anterior panel provides additional abdominal lightweight comfort and support!
- All panels are easily and quickly removed with no tools.



Posterior Panel



Velcro Strap

Anterior Panel Attachment

S.P.I.N.E. Plus

Clinical Guidelines:

"CYBERTECH Medical braces are medical products registered with the FDA as a "medical device". CYBERTECH Medical does not prescribe braces or determine the appropriateness of any particular bracing methodology for the individual. Medical advice should be sought on the appropriateness and the condition under which a given brace should be worn".

Fitting Specifications:

Size	Circ.	Posterior Pnl.	Anterior Pnl.	Item#: SPINE Plus	Item#: SPINE Sport
XSmall	20"-25"	7 1/2"	4 3/4"	ESFP-NB-0	
Small	25"- 30"	7 1/2"	4 3/4"	ESFP-NB-1	OM500-S
Medium	30"-35"	8"	5 1/2"	ESFP-NB-2	OM500-M
Large	35"-40"	8 1/2"	6"	ESFP-NB-3	OM500-L
X-Large	40"-45"	9"	6 7/8"	ESFP-NB-4	OM500-XL
XX-Large	45"-50"	9 1/2"	7 3/4"	ESFP-NB-5	OM500-XXL
XXX-Large	50"-55"	10"	9"	ESFP-NB-6	OM500-XXXL

Braces come in Black color. Please use "NW" for white color. XS size comes in Black only.

U.S. Patent No. RE35,940

Brace come in Black color.



Fitting Instructions:

- Fasten the handle near the pulley unit. The handle is always on the left.
- Stand holding the soft side panel apart, centering the back panel on the spine.
- Wrap the left side to the front positioning the front panel.
- Wrap the right side on top of the left panel to close and attach in front. Adjust the front or back panels as needed.
- Pull the handle out and away from the power unit to desired compression and support. Place the handle on the fabric fastener band securely.

Tips:

- Be sure to start the handle near the pulley unit.
- Extend or pull out the brace before application.
- Adjust tension for sitting or standing, experiment to meet your needs.
- Hand wash in warm water, rinse well, towel and air dry.

Posterior and Anterior Panels:

- Dynamic Posterior panels conform to each patients individual body structure for improved comfort and support
- Anterior panels are easily attached with integrated Velcro strap.

Dynamic Posterior panel easy and secure attachment and removal:



1. Lift up attachment strap



2. Place standard rivet head in slot



3. Snap into place in slot



Placed rivet head



4. Slide rivet head to end of slot to remove



1815 Wright Ave., La Verne, CA 91750
Ph: 1.800.220.4224 Fax: 1.909.447.7058
www.cybertechmedical.com

© 2011 Bio Cybernetics International, Inc. All rights reserved. Rev. 03/11



Authorized Rep.
RSQA / RSQR Ltd.
www.rsqa.co.uk
info@rsqa.co.uk