

Dr. Zen Therapeutic Shoes Product Information and Sizing Chart

2/3/2011

TWS - Triple Width System

Medium is B, C, D Wide is E, 2E, 3E X-Wide is 4E, 5E, 6E Each Triple Width System (**TWS**) shoe has 3 removable inserts, every time you remove an insert the shoe becomes 1 width wider. The platforms are designed specifically for the widest width. Every added insert takes up one width. This gives the fitter, patient, etc the option to customize the fit of the shoe to accommodate swelling, pronation, certain deformities where one foot is wider or different in some way as opposed to the other foot, or where one foot is longer or shorter then the other. Triple Width System shoes are indicated with a **TWS** next to their name below.

Aries - TWS			• TWS - Triple Width System • Lace & Velcro Closure • Widths available from B - 6E • Reinforced Counter • Achilles Tendon Support • Breathable Leather / Mesh • Thomas Heel • 3/8 Extra Depth • Seamless Construction • Rocker Sole
------------------------	---	--	--

Product #	Shoe	Color	Closure	Width	Women's Size	Men's Size
745-01	Aries	Grey	Velcro	Medium (B, C, D)	8-13,14,15,16	7-12,13,14,15
745-02	Leather / Mesh			Wide (E, 2E, 3E)	8-13,14,15,16	7-12,13,14,15
745-03				X-Wide (4E, 5E, 6E)	8-13,14,15,16	7-12,13,14,15
755-01	Aries	Grey	Lace	Medium (B, C, D)	8-13,14,15,16	7-12,13,14,15
755-02	Leather / Mesh			Wide (E, 2E, 3E)	8-13,14,15,16	7-12,13,14,15
755-03				X-Wide (4E, 5E, 6E)	8-13,14,15,16	7-12,13,14,15

Gemini - TWS			• TWS - Triple Width System • Lace & Velcro Closure • Widths available from B - 6E • Reinforced Counter • Achilles Tendon Support • Breathable Leather / Mesh • Thomas Heel • 3/8 Extra Depth • Seamless Construction • Rocker Sole
-------------------------	---	--	--

Product #	Shoe	Color	Closure	Width	Women's Size	Men's Size
645-01	Gemini	Grey	Velcro	Medium (B, C, D)	8-13,14,15,16	7-12,13,14,15
645-02	Leather / Mesh			Wide (E, 2E, 3E)	8-13,14,15,16	7-12,13,14,15
645-03				X-Wide (4E, 5E, 6E)	8-13,14,15,16	7-12,13,14,15
655-01	Gemini	Grey	Lace	Medium (B, C, D)	8-13,14,15,16	7-12,13,14,15
655-02	Leather / Mesh			Wide (E, 2E, 3E)	8-13,14,15,16	7-12,13,14,15
655-03				X-Wide (4E, 5E, 6E)	8-13,14,15,16	7-12,13,14,15

Dr. Zen Therapeutic Shoes
Product Information and Sizing Chart

2/3/2011

Storm - TWS			• TWS - Triple Width System • Lace & Velcro Closure • Widths available from B - 6E • Reinforced Counter • Achilles Tendon Support • Breathable Leather / Mesh • Thomas Heel • 3/8 Extra Depth • Seamless Construction • Rocker Sole
------------------------------	---	--	--

Product #	Shoe	Color	Closure	Width	Women's Size	Men's Size
545-01	Storm	Black / Blue	Velcro	Medium (B, C, D)	8-13,14,15,16	7-12,13,14,15
545-02	Leather / Mesh			Wide (E, 2E, 3E)	8-13,14,15,16	7-12,13,14,15
545-03				X-Wide (4E, 5E, 6E)	8-13,14,15,16	7-12,13,14,15
555-01	Storm	Black / Blue	Lace	Medium (B, C, D)	8-13,14,15,16	7-12,13,14,15
555-02	Leather / Mesh			Wide (E, 2E, 3E)	8-13,14,15,16	7-12,13,14,15
555-03				X-Wide (4E, 5E, 6E)	8-13,14,15,16	7-12,13,14,15

Lori - TWS	 	 	• TWS - Triple Width System • Lace & Velcro Closure • Widths available from B - 6E • Reinforced Counter • Achilles Tendon Support • Breathable Leather / Mesh • Thomas Heel • 3/8 Extra Depth • Seamless Construction • Rocker Sole
-----------------------------	--	--	--

Product #	Shoe	Color	Closure	Width	Women's Size	Men's Size
945-01	Lori	White / Pink	Velcro	Medium (B, C, D)	6-11,12	5-10,11
945-02	Leather / Mesh			Wide (E, 2E, 3E)	6-11,12	5-10,11
945-03				X-Wide (4E, 5E, 6E)	6-11,12	5-10,11
955-01	Lori	White / Pink	Lace	Medium (B, C, D)	6-11,12	5-10,11
955-02	Leather / Mesh			Wide (E, 2E, 3E)	6-11,12	5-10,11
955-03				X-Wide (4E, 5E, 6E)	6-11,12	5-10,11
946-01	Lori	White / Lavendar	Velcro	Medium (B, C, D)	6-11,12	5-10,11
946-02	Leather / Mesh			Wide (E, 2E, 3E)	6-11,12	5-10,11
946-03				X-Wide (4E, 5E, 6E)	6-11,12	5-10,11
956-01	Lori	White / Lavendar	Lace	Medium (B, C, D)	6-11,12	5-10,11
956-02	Leather / Mesh			Wide (E, 2E, 3E)	6-11,12	5-10,11
956-03				X-Wide (4E, 5E, 6E)	6-11,12	5-10,11

Dr. Zen Therapeutic Shoes
Product Information and Sizing Chart

2/3/2011

**Anne -
TWS**



- TWS - Triple Width System
- Velcro Closure
- Widths available from B - 6E
- Reinforced Counter
- Achilles Tendon Support
- Breathable Leather / Mesh
- Thomas Heel
- 3/8 Extra Depth
- Seamless Construction
- Rocker Sole

Product #	Shoe	Color	Closure	Width	Women's Size	Men's Size
222-01	Anne	Black	Velcro	Medium (B, C, D)	6-11,12	5-10,11
222-02	Leather			Wide (E, 2E, 3E)	6-11,12	5-10,11
222-03				X-Wide (4E, 5E, 6E)	6-11,12	5-10,11
221-01	Anne	Brown	Velcro	Medium (B, C, D)	6-11,12	5-10,11
221-02	Leather			Wide (E, 2E, 3E)	6-11,12	5-10,11
221-03				X-Wide (4E, 5E, 6E)	6-11,12	5-10,11

**Beth -
TWS**



- TWS - Triple Width System
- Velcro Closure
- Widths available from B - 6E
- Reinforced Counter
- Achilles Tendon Support
- Breathable Leather / Mesh
- Thomas Heel
- 3/8 Extra Depth
- Seamless Construction
- Rocker Sole

Product #	Shoe	Color	Closure	Width	Women's Size	Men's Size
178-01	Beth	Black	Velcro	Medium (B, C, D)	6-11,12	5-10,11
178-02	Leather			Wide (E, 2E, 3E)	6-11,12	5-10,11
178-03				X-Wide (4E, 5E, 6E)	6-11,12	5-10,11

**Chris -
TWS**



- TWS - Triple Width System
- Velcro Closure
- Widths available from B - 6E
- Reinforced Counter
- Achilles Tendon Support
- Breathable Leather / Mesh
- Thomas Heel
- 3/8 Extra Depth
- Seamless Construction
- Rocker Sole

Product #	Shoe	Color	Closure	Width	Women's Size	Men's Size
799-01	Chris	Brown	Velcro	Medium (B, C, D)	6-11,12	5-10,11
799-02	Leather			Wide (E, 2E, 3E)	6-11,12	5-10,11
799-03				X-Wide (4E, 5E, 6E)	6-11,12	5-10,11

Dr. Zen Therapeutic Shoes
Product Information and Sizing Chart

2/3/2011

**Linda -
TWS**



- TWS - Triple Width System
- Velcro Closure
- Widths available from B - 6E
- Reinforced Counter
- Achilles Tendon Support
- Breathable Leather / Mesh
- Thomas Heel
- 3/8 Extra Depth
- Seamless Construction
- Rocker Sole

Product #	Shoe	Color	Closure	Width	Women's Size	Men's Size
444-01	Linda	Brown	Velcro	Medium (B, C, D)	6-11,12	5-10,11
444-02	Leather			Wide (E, 2E, 3E)	6-11,12	5-10,11
444-03				X-Wide (4E, 5E, 6E)	6-11,12	5-10,11

**Melissa -
TWS**



- TWS - Triple Width System
- Velcro Closure
- Widths available from B - 6E
- Reinforced Counter
- Achilles Tendon Support
- Breathable Leather / Mesh
- Thomas Heel
- 3/8 Extra Depth
- Seamless Construction
- Rocker Sole

Product #	Shoe	Color	Closure	Width	Women's Size	Men's Size
111-01	Melissa	Black	Velcro	Medium (B, C, D)	6-11,12	5-10,11
111-02	Mesh			Wide (E, 2E, 3E)	6-11,12	5-10,11
111-03				X-Wide (4E, 5E, 6E)	6-11,12	5-10,11

**Nellie -
TWS**



- TWS - Triple Width System
- Lace Closure
- Widths available from B - 6E
- Reinforced Counter
- Achilles Tendon Support
- Breathable Leather / Mesh
- Thomas Heel
- 3/8 Extra Depth
- Seamless Construction
- Rocker Sole

Product #	Shoe	Color	Closure	Width	Women's Size	Men's Size
402-01	Nellie	Black	Lace	Medium (B, C, D)	6-11,12	5-10,11
402-02	Leather			Wide (E, 2E, 3E)	6-11,12	5-10,11
402-03				X-Wide (4E, 5E, 6E)	6-11,12	5-10,11

Dr. Zen Therapeutic Shoes
Product Information and Sizing Chart

2/3/2011

Stella - TWS			• TWS - Triple Width System • Velcro Closure • Widths available from B - 6E • Reinforced Counter • Achilles Tendon Support • Breathable Leather / Mesh • Thomas Heel • 3/8 Extra Depth • Seamless Construction • Rocker Sole
-------------------------------	---	--	---

<i>Product #</i>	<i>Shoe</i>	<i>Color</i>	<i>Closure</i>	<i>Width</i>	<i>Women's Size</i>	<i>Men's Size</i>
855-01	Stella	Black	Velcro	Medium (B, C, D)	6-11,12	5-10,11
855-02	Leather / Mesh			Wide (E, 2E, 3E)	6-11,12	5-10,11
855-03				X-Wide (4E, 5E, 6E)	6-11,12	5-10,11

Valerie - TWS			• TWS - Triple Width System • Velcro Closure • Widths available from B - 6E • Reinforced Counter • Achilles Tendon Support • Breathable Leather / Mesh • Thomas Heel • 3/8 Extra Depth • Seamless Construction • Rocker Sole
--------------------------------	--	---	---

<i>Product #</i>	<i>Shoe</i>	<i>Color</i>	<i>Closure</i>	<i>Width</i>	<i>Women's Size</i>	<i>Men's Size</i>
125-01	Valerie	Black	Velcro	Medium (B, C, D)	6-11,12	5-10,11
125-02	Mesh / Canvas			Wide (E, 2E, 3E)	6-11,12	5-10,11
125-03				X-Wide (4E, 5E, 6E)	6-11,12	5-10,11
126-01	Valerie	Silver	Velcro	Medium (B, C, D)	6-11,12	5-10,11
126-02	Mesh / Canvas			Wide (E, 2E, 3E)	6-11,12	5-10,11
126-03				X-Wide (4E, 5E, 6E)	6-11,12	5-10,11

Dr. Zen Therapeutic Shoes
Product Information and Sizing Chart

2/3/2011

Ben Boot - TWS			• TWS - Triple Width System • Lace Closure • Widths available from B - 6E • Reinforced Counter • Achilles Tendon Support • Breathable Leather / Mesh • Thomas Heel • 3/8 Extra Depth • Seamless Construction • Rocker Sole
---------------------------	---	--	---

Product #	Shoe	Color	Closure	Width	Women's Size	Men's Size
900-01	Ben	Black	Lace	Medium (B, C, D)	8-13,14,15,16	7-12,13,14,15
900-02	Leather			Wide (E, 2E, 3E)	8-13,14,15,16	7-12,13,14,15
900-03				X-Wide (4E, 5E, 6E)	8-13,14,15,16	7-12,13,14,15

Charles - TWS			• TWS - Triple Width System • Lace Closure • Widths available from B - 6E • Reinforced Counter • Achilles Tendon Support • Breathable Leather / Mesh • Thomas Heel • 3/8 Extra Depth • Seamless Construction • Rocker Sole
--------------------------	--	---	---

Product #	Shoe	Color	Closure	Width	Women's Size	Men's Size
331-01	Charles	Black	Lace	Medium (B, C, D)	8-13,14,15,16	7-12,13,14,15
331-02	Leather			Wide (E, 2E, 3E)	8-13,14,15,16	7-12,13,14,15
331-03				X-Wide (4E, 5E, 6E)	8-13,14,15,16	7-12,13,14,15
332-01	Charles	Brown	Lace	Medium (B, C, D)	8-13,14,15,16	7-12,13,14,15
332-02	Leather			Wide (E, 2E, 3E)	8-13,14,15,16	7-12,13,14,15
332-03				X-Wide (4E, 5E, 6E)	8-13,14,15,16	7-12,13,14,15
333-01	Charles	Black	Velcro	Medium (B, C, D)	8-13,14,15,16	7-12,13,14,15
333-02	Leather			Wide (E, 2E, 3E)	8-13,14,15,16	7-12,13,14,15
333-03				X-Wide (4E, 5E, 6E)	8-13,14,15,16	7-12,13,14,15
334-01	Charles	Brown	Velcro	Medium (B, C, D)	8-13,14,15,16	7-12,13,14,15
334-02	Leather			Wide (E, 2E, 3E)	8-13,14,15,16	7-12,13,14,15
334-03				X-Wide (4E, 5E, 6E)	8-13,14,15,16	7-12,13,14,15

Dr. Zen Therapeutic Shoes
Product Information and Sizing Chart

2/3/2011

Legend Golf - TWS			• Golf shoe • Creme Lace • Extra depth for added protection • Rubber sole with cleats
--------------------------	---	--	--

Product #	Shoe	Color	Closure	Width	Women's Size	Men's Size
418-01	Legend Golf	Crème	Lace	Medium (B, C, D)	7-13	6-12
418-03	Leather			X-Wide (D, E, 2E)	7-13	6-12
418-05				XXX-Wide (2E, 3E, 4E)	7-13,14	6-12,13

Legend - TWS			• Leather dress and walking shoe • Lace closure • Rubber sole
---------------------	--	---	---

Product #	Shoe	Color	Closure	Width	Women's Size	Men's Size
310-01	Legend	Black	Lace	Medium (B, C, D)	7-13,14	6-12,13
310-03	Leather			X-Wide (D, E, 2E)	7-13,14	6-12,13
310-05				XXX-Wide (2E, 3E, 4E)	7-13,14	6-12,13
311-00	Legend	Crème	Lace	Narrow (AA, A, B)	7-11	6-10
311-01	Leather			Medium (B, C, D)	7-13	6-12
311-03				X-Wide (D, E, 2E)	7-13,14	6-12,13

Nautica Boat - TWS			• Leather boat shoe • Lace closure • Extra depth for added protection • Rubber sole
---------------------------	---	--	--

Product #	Shoe	Color	Closure	Width	Women's Size	Men's Size
615-01	Nautica	Brown	Lace	Medium (B, C, D)	7-12,13	6-11,12
615-03	Leather			X-Wide (D, E, 2E)	7-12,13,14	6-11,12,13
615-05				XXX-Wide (2E, 3E, 4E)	8-12,13,14	7-11,12,13

Dr. Zen Therapeutic Shoes
Product Information and Sizing Chart

2/3/2011

Dawn			• Soft Leather • Black with Velcro closure • Narrow, Med, Wide & X-Wide • Rubber sole • Kosmo has similar shoe pattern
-------------	---	--	--

Product #	Shoe	Color	Closure	Width	Women's Size	Men's Size
248-00	Dawn	Black	Velcro	Narrow (A)	5-11	4-10
248-01	Leather			Medium (A-D)	5-12	4-11
248-02				Wide (E)	5-12	4-11
248-03				X-Wide (EE)	5-12	4-11
249-00	Dawn	Beige	Velcro	Narrow (A)	5-10	4-9
249-01	Leather			Medium (A-D)	5-11	4-10
249-02				Wide (E)	5-12	4-11
249-03				X-Wide (EE)	6-12	5-11

Kosmo			• Soft Leather • Black with Lace closure • Wide, X-Wide, XX-Wide, 4E and 6E • Rubber sole • Dawn has similar shoe pattern
--------------	--	---	---

Product #	Shoe	Color	Closure	Width	Women's Size	Men's Size
899-02	Kosmo	Black	Lace	Wide (E)	8-12,13,14	7-11,12,13
899-03	Leather			X-Wide (EE)	8-12,13,14,15	7-11,12,13,14
899-04				XX-Wide (3E)	9,10,11,12,13,14,15	8,9,10,11,12,13,14
899-05				XXX-Wide (4E)	8, 9, 10,11,12,13,14,15	7, 8, 9,10,11,12,13,14
899-06				XXXXX-Wide (6E)	8, 9, 10,11,12,13,14,15	7, 8, 9,10,11,12,13,14

Dr. Zen Therapeutic Shoes
Product Information and Sizing Chart

2/3/2011

Gator			• Leather • Black, Brown, Beige • Velcro or Lace closure • Med, Wide and X-Wide • Rubber sole • Top stitching creates an inside seam
--------------	---	--	---

<i>Product #</i>	<i>Shoe</i>	<i>Color</i>	<i>Closure</i>	<i>Width</i>	<i>Women's Size</i>	<i>Men's Size</i>
195-01	Gator	Black	Velcro	Medium (A-D)	7-12,13	6-11,12
195-02	Leather			Wide (E)	6-12,13,14,15	5-11,12,13,14
195-03				X-Wide (EE)	8-12,13,14,15	7-11,12,13,14
196-01	Gator	Black	Lace	Medium (A-D)	7-12,13	6-11,12
196-02	Leather			Wide (E)	6-12,13,14,15	5-11,12,13,14
196-03				X-Wide (EE)	8-12,13,14,15	7-11,12,13,14
295-01	Gator	Brown	Velcro	Medium (A-D)	7-12,13	6-11,12
295-02	Leather			Wide (E)	6-12,13,14,15	5-11,12,13,14
295-03				X-Wide (EE)	8-12,13,14,15	7-11,12,13,14
296-01	Gator	Brown	Lace	Medium (A-D)	7-12,13	6-11,12
296-02	Leather			Wide (E)	6-12,13,14,15	5-11,12,13,14
296-03				X-Wide (EE)	8-12,13,14,15	7-11,12,13,14
395-01	Gator	Beige	Velcro	Medium (A-D)	7-12,13	6-11,12
395-02	Leather			Wide (E)	6-12,13,14,15	5-11,12,13,14
395-03				X-Wide (EE)	8-12,13,14,15	7-11,12,13,14
396-01	Gator	Beige	Lace	Medium (A-D)	7-12,13	6-11,12
396-02	Leather			Wide (E)	6-12,13,14,15	5-11,12,13,14
396-03				X-Wide (EE)	8-12,13,14,15	7-11,12,13,14

Dr. Zen Therapeutic Shoes
Product Information and Sizing Chart

2/3/2011

Jenny			• Black with Leather toe area & Velcro closure • Narrow, Med, Wide and X-Wide • Taupe with Zennon toe area & Velcro closure • Rubber sole • Shoe design helps alleviate top swelling • Zennon accommodates swelling & toe deformities
--------------	---	--	--

Product #	Shoe	Color	Closure	Width	Women's Size	Men's Size
148-00	Jenny	Black	Velcro	Narrow (A)	5-11	4-10
148-01	Leather			Medium (A-D)	6-10	5-9
148-02				Wide (E)	5-11,12	4-10,11
148-03				X-Wide (EE)	7-11,12	6-10,11
158-00	Jenny	Taupe	Velcro	Narrow (A)	5-11	4-10
158-01	Zennon			Medium (A-D)	6-10	5-9
158-02				Wide (E)	5-11,12	4-10,11
158-03				X-Wide (EE)	7-11,12	6-10,11
168-00	Jenny	Black	Velcro	Narrow (A)	5-11	4-10
168-01	Zennon			Medium (A-D)	6-10	5-9
168-02				Wide (E)	5-11,12	4-10,11
168-03				X-Wide (EE)	7-11,12	6-10,11

Leather			• Velcro & Lace* closure • Large toe box for deformities • Med, Wide, X-Wide, XX-Wide, 4E and 6E • Rubber sole with wide heel and rocker toe • High volume / space throughout shoe • High - thick toe box for added protection • Zennon has similar shoe pattern * Lace is a Velcro/Lace closure
----------------	---	--	---

Product #	Shoe	Color	Closure	Width	Women's Size	Men's Size
400-01	Leather	Black	Velcro	Medium (A-D)	4-12,13,14	3-11,12,13
400-02				Wide (E)	4-12,13,14	3-11,12,13
400-03				X-Wide (EE)	5-12,13,14,15	4-11,12,13,14
400-04				XX-Wide (3E)	8,9,10,11,12,13,14,15,16	7,8,9,10,11,12,13,14,15
400-05				XXX-Wide (4E)	8,9,10,11,12,13,14,15	7,8,9,10,11,12,13,14
400-06				XXXXX-Wide (6E)	8,9,10,11,12,13,14,15	7,8,9,10,11,12,13,14
401-01	Leather	Black	Lace	Medium (A-D)	4-12,13,14	3-11,12,13
401-02				Wide (E)	4-12,13,14,15	3-11,12,13,14
401-03				X-Wide (EE)	5-12,13,14,15	4-11,12,13,14
401-04				XX-Wide (3E)	8,9,10,11,12,13,14,15,16	7,8,9,10,11,12,13,14,15
401-05				XXX-Wide (4E)	8,9,10,11,12,13,14,15	7,8,9,10,11,12,13,14
401-06				XXXXX-Wide (6E)	8,9,10,11,12,13,14,15	7,8,9,10,11,12,13,14

Dr. Zen Therapeutic Shoes
Product Information and Sizing Chart

2/3/2011

New Yorker			• Leather dress shoe with Wing-Tip stitching • Black Lace closure • Wide, X-Wide, XX-Wide, 4E and 6E • High - thick toe box for added protection • Rubber sole
-------------------	---	--	--

Product #	Shoe	Color	Closure	Width	Women's Size	Men's Size
999-02	New Yorker	Black	Lace	Wide (E)	8-12,13,14,15	7-11,12,13,14
999-03	Leather			X-Wide (EE)	8-12,13,14,15	7-11,12,13,14
999-04				XX-Wide (3E)	8,9,10,11,12,13,14,15	7,8,9,10,11,12,13,14
999-05				XXX-Wide (4E)	8,9,10,11,12,13,14,15	7,8,9,10,11,12,13,14
999-06				XXXXX-Wide (6E)	8,9,10,11,12,13,14,15	7,8,9,10,11,12,13,14

Raven			• Leather with partial Zennon top • Velcro closure • Med, Wide and X-Wide • Rubber sole • Top stitching creates an inside seam
--------------	--	---	--

Product #	Shoe	Color	Closure	Width	Women's Size	Men's Size
114-01	Raven	Beige	Velcro	Medium (A-D)	6-12	5-11
114-02	Leather			Wide (E)	5-12,13,14	4-11,12,13
114-03				X-Wide (EE)	7-12,13,14	6-11,12,13
113-01	Raven	Black	Velcro	Medium (A-D)	6-12	5-11
113-02	Leather			Wide (E)	5-12,13,14,15	4-11,12,13,14
113-03				X-Wide (EE)	7-12,13,14,15	6-11,12,13,14

Sandal		NOT a Medicare Covered Item - No Inserts are available		• Not covered by Medicare • Velcro closure • Medium only • No Inserts available • Light Brown - Women • Dark Brown - Mens
---------------	---	--	--	--

Product #	Shoe	Color	Closure	Width	Women's Size	Men's Size
712-02	Sandal	Tan	Velcro	Wide (E)	7,8,9,10,11,12,13,14	6,7,8,9,10,11,12,13
711-02		Dark Brown	Velcro	Wide (E)	6,7,8,9,10,11,12	5,6,7,8,9,10,11

Dr. Zen Therapeutic Shoes
Product Information and Sizing Chart

2/3/2011

Sport I			• Leather • White or Black • Velcro or Lace • Med, Wide and X-Wide • Non-Slip Rubber sole • Exercise / Walking shoe • High back for narrow heel patients • Padded tongue and heel area • Sport II has similar shoe pattern
----------------	---	--	--

Product #	Shoe	Color	Closure	Width	Women's Size	Men's Size
198-01	Sport I	White	Velcro	Med (A-D)	7-10	6-9
198-02	Leather			Wide (E)	5-12,13,14,15	4-11,12,13,14
198-03				X-Wide (EE)	6.5-12,13,14	5.5-11,12,13
199-01	Sport I	White	Lace	Med (A-D)	7-10	6-9
199-02	Leather			Wide (E)	5-12,13,14,15	4-11,12,13,14
199-03				X-Wide (EE)	6.5-12,13,14	5.5-11,12,13
298-01	Sport I	Black	Velcro	Med (A-D)	7-10	6-9
298-02	Leather			Wide (E)	5-12,13,14,15	4-11,12,13,14
298-03				X-Wide (EE)	6.5-12,13,14	5.5-11,12,13
299-01	Sport I	Black	Lace	Med (A-D)	7-10	6-9
299-02	Leather			Wide (E)	5-12,13,14,15	4-11,12,13,14
299-03				X-Wide (EE)	6.5-12,13,14	5.5-11,12,13

Sport II			• Leather • White or Black • Velcro or Lace • Med, Wide and X-Wide • Non-Slip Rubber sole • Exercise / Walking shoe • High back for narrow heel patients • Padded tongue and heel area • Sport I has similar shoe pattern
-----------------	---	--	---

Product #	Shoe	Color	Closure	Width	Women's Size	Men's Size
398-01	Sport II	White	Velcro	Medium (A-D)	5-12,13,14,15	4-11,12,13,14
398-02	Leather			Wide (E)	5-12,13,14,15	4-11,12,13,14
398-03				X-Wide (EE)	6.5-12,13,14,15	5.5-11,12,13,14
399-01	Sport II	White	Lace	Medium (A-D)	5-12,13	4-11,12
399-02	Leather			Wide (E)	5-12,13,14,15	4-11,12,13,14
399-03				X-Wide (EE)	6.5-12,13,14,15	5.5-11,12,13,14
498-01	Sport II	Black	Velcro	Medium (A-D)	5-12,13	4-11,12
498-02	Leather			Wide (E)	5-12,13,14,15	4-11,12,13,14
498-03				X-Wide (EE)	6.5-12,13,14,15	5.5-11,12,13,14
499-01	Sport II	Black	Lace	Medium (A-D)	5-12,13	4-11,12
499-02	Leather			Wide (E)	5-12,13,14,15	4-11,12,13,14
499-03				X-Wide (EE)	6.5-12,13,14,15	5.5-11,12,13,14

Dr. Zen Therapeutic Shoes
Product Information and Sizing Chart

2/3/2011

Suede			• Leather & Suede • Brown suede with Lace / Velcro closure • Wide • Rubber sole with wide heel and rocker toe • High volume / space throughout shoe • High - thick toe box for added protection • Zennon & Leather have similar shoe pattern
--------------	---	--	--

Product #	Shoe	Color	Closure	Width	Women's Size	Men's Size
801-01	Suede	Brown	Lace	Medium (A-D)	6.5-8.5	5.5-7.5
801-02				Wide (E)	5-12,13,14	4-11,12,13
801-03				X-Wide (EE)	8-10	7-9

Titan			• Fabric & Zennon material • Med, Wide, and X-Wide • Gently wraps the foot like a bedroom slipper • Indoor shoe only • Washable, bag included • Zennon covering alleviates toe deformities
--------------	---	--	---

Product #	Shoe	Color	Closure	Width	Women's Size	Men's Size
214-01	Titan	Beige	Velcro	Medium (A-D)	6-15	5-14
214-02	Lycra			Wide (E)	6-15	5-14
214-03				X-Wide (EE)	6-15,16,17	5-14,15,16
213-01	Titan	Black	Velcro	Medium (A-D)	5-15	4-14
213-02	Lycra			Wide (E)	5-15	4-14
213-03				X-Wide (EE)	5-15,16,17	4-14,15,16
213-06				XXXXX-Wide (6E)	6-14,15	5-13, 14

Work Boot			• Leather • Black with Lace / Velcro • Wide, X-Wide, XX-Wide • Hi-top provides ankle support • Rubber sole with wide heel and rocker toe • High volume / space throughout shoe • High - thick toe box for added protection • Zennon & Leather have similar shoe pattern
------------------	---	--	--

Product #	Shoe	Color	Closure	Width	Women's Size	Men's Size
901-02	Work Boot	Black	Lace/Velcro	Wide (E)	7-12,13,14	6-11,12,13
901-03	Leather			X-Wide (EE)	9-12,13,14,15	8-11,12,13,14
901-04				XX-Wide (3E)	9,10,11,12,13,14,15	8,9,10,11,12,13,14

Dr. Zen Therapeutic Shoes
Product Information and Sizing Chart

2/3/2011

Zennon			• Leather with Zennon toe area • Velcro or Lace • Med, Wide, X-Wide, XX-Wide, 4E and 6E • Rubber sole with wide heel and rocker toe • High volume / space throughout shoe • High - thick toe box for added protection • Leather has similar shoe pattern • Zennon accommodates swelling & toe deformities
---------------	---	--	--

Product #	Shoe	Color	Closure	Width	Women's Size	Men's Size
300-01	Zennon	Black	Velcro	Medium (A-D)	4-12,13,14	3-11,12,13
100-01	Lycra			Wide (E)	4-12,13,14	3-11,12,13
200-01				X-Wide (EE)	4-12,13,14,15	3-11,12,13,14
100-04				XX-Wide (3E)	8,9,10,11,12,13,14,15,16	7,8,9,10,11,12,13,14,15
100-05				XXX-Wide (4E)	8,9,10,11,12,13,14,15	7,8,9,10,11,12,13,14
100-06				XXXXX-Wide (6E)	8,9,10,11,12,13,14,15	7,8,9,10,11,12,13,14
301-01	Zennon	Black	Lace	Medium (A-D)	5-10	4-9
101-01	Lycra			Wide (E)	4-12,13,14	3-11,12,13
201-01				X-Wide (EE)	4-12,13,14,15,16	3-11,12,13,14,15
101-04				XX-Wide (3E)	8,9,10,11,12,13,14,15,16	7,8,9,10,11,12,13,14,15
101-05				XXX-Wide (4E)	8,9,10,11,12,13,14,15	7,8,9,10,11,12,13,14
101-06				XXXXX-Wide (6E)	8,9,10,11,12,13,14,15	7,8,9,10,11,12,13,14
300-03	Zennon	Taupe	Velcro	Medium (A-D)	4-12,13,14	3-11,12,13
100-03	Lycra			Wide (E)	4-12,13,14	3-11,12,13
200-03				X-Wide (EE)	4-12,13,14	3-11,12,13
300-02	Zennon	Brown	Velcro	Medium (A-D)	5-10	4-9
100-02	Lycra			Wide (E)	4-12,13,14	3-11,12,13
200-02				X-Wide (EE)	8-12,13,14	7-11,12,13

Dr. Zen Therapeutic Shoes
Product Information and Sizing Chart

2/3/2011

Gauntlet		Gauntlet Ankle Foot Orthotic - AFO ▪ Fabric with Velcro straps ▪ Covers toes, open at heel ▪ Ankle & Joint Support ▪ Defective Circulation ▪ Joint Stiffness ▪ Joint Swelling ▪ Arthritis ▪ Disuse Atrophy ▪ Keeps feet warm at night
-----------------	---	--

	Product			Size	Women's Size	Men's Size
	Gauntlet					
				X - Small	4 - 5	3 - 4
				Small	6 - 7	5 - 6
				Medium	8 - 9	7 - 8
				Large	10 - 11	9 - 10
				X-Large	12 & Up	11 & Up

Diabetic Socks		Diabetic Socks ▪ Lightweight ▪ Soft cotton 88% cotton 12% spandex ▪ Seamless socks ▪ Not reimbursable by Medicare
-----------------------	--	---

	Product			Size	Women's Size	Men's Size
	Socks					
				Small	4 - 6	3 - 5
				Medium	7 - 9	6 - 8
				Large	10 - 11	9 - 10
				X-Large	12 & Up	11 & Up

Dr. Zen Therapeutic Shoes
Product Information and Sizing Chart

2/3/2011

**Custom Therapeutic
Shoe Inserts**



Custom 3-Layer Inserts Materials: Top layer is fine celled Pro-Cell Eva, the insert will naturally heat-mold to the foot without abrasion. The center layer is Pron Foam that will not bottom out and continue to provide a comfortable cushion. The bottom layer is a firmer Pro-cell Plus Eva designed to hold the shape of the insole.

Specifications: Designed specifically for patient's foot and selected shoe. Insert works together with shoe to support the foot and arch, as well as protecting the foot's sole. Insert is specifically designed to create less friction on the foot and fits inside the shoes.

**Impression Box &
Ritz Stick**



Impression box is used to obtain a 3-D impression of the patient's feet. The custom inserts are created from this impression. Obtaining a deep impression of the patient's feet along with accurate measurements is very important to create an accurate custom insert.

Ritz stick is used to obtain adequate foot measurements and is vital to fitting shoes properly.