A revolutionary way to reduce pain and swelling through the hot, cold, contrast and compression therapies you trust.

put a wrap on pain



# Physicians and patients are fighting pain with precision, ease and flexibility.

Real relief, right where you need it. Reduce pain and swelling anywhere in the body-without drugs-in one breakthrough device. The **SADMERC** - approved VitalWrap System with its innovative family of wraps delivers targeted relief through hot, cold, contrast and compression therapies with precision and comfort.

#### SIMPLE:

Portable, easy-to-use design for patient self-treatment.

**PRECISE:** Control temperature with one dial, 40° - 105°F.

**FLEXIBLE:** Hot, cold, contrast and compression therapies.

COMPLETE:

Family of therapeutic wraps treats multiple body areas.

**Wide-ranging applications — One versatile solution** Arthritis • Back pain • Sprains, strains and inflammation • Post-surgical recovery Pain management • Functional mobility • Many other conditions





put a wrap on pain

# Systems order no.

00002	2 VitalWrap System 120 V <sub>ac</sub> , 60 Hz (control unit, tubing set, standard VW)	
00013	VitalWrap System 120 V <sub>aC</sub> , 60 Hz (control unit, tubing set, without standard VW)	
Wraps	5	
00255	Medium VitalWrap (4.25" x 52")	for elbow and calf
00001	Standard VitalWrap (6" x 60")	general purpose wrap for calf, knee, hamstring, hip, back, chest and abdomen
00118	Large VitalWrap (8" x 72")	broad coverage wrap for leg and back
00388	Back VitalWrap	for use on back
00119	Medium Shoulder VitalWrap	for use on either shoulder
00395	Large Shoulder VitalWrap	for use on either shoulder
00394	Wrist VitalWrap	for hand and wrist
00396	Hand VitalWrap	for arthritic hands
00397	Medium Knee VitalWrap	for post-surgery and CPM applications
00392	Full Knee VitalWrap	for post-surgery and CPM applications
00398	Ankle VitalWrap	for foot and ankle
00122	Cervical VitalWrap	for neck applications
00967	Hip VitalWrap	for use on hip
01061	Cap VitalWrap	for use on head

#### VitalPad<sup>®</sup> made exclusively for XBACK<sup>®</sup> Solution LSO Brace

01311	VitalPad for XBACK Solution, S	for XBACK Solution Model T102	
01315	VitalPad for XBACK Solution, M	for XBACK Solution Model T102	
01319	VitalPad for XBACK Solution, L	for XBACK Solution Model T102	
01323	VitalPad for XBACK Solution, XL	for XBACK Solution Model T102	
01327	VitalPad for XBACK Solution, 2XL	for XBACK Solution Model T102	
01331	VitalPad for XBACK Solution, 3XL	for XBACK Solution Model T102	
01335	VitalPad for XBACK Solution, 4XL	for XBACK Solution Model T102	
XBACK <sup>®</sup> Bracing Services, Inc. is registered trademark.			

#### Accessories

00094 Carrying Case	for easy portability of the VW System
00034 VitalWrap Tubing Set	
00987 VitalWrap "Y" Tubing Set	for use on 2 limbs
00262 Hospital Grade Power Cord	

### SADMERC Durable Medical Equipment Codes code

E0217	water circulating heat pad with	pump
E0218	Water circulating cold pad with	pump

E0249 Pad for water circulation heat unit

VitalWear, Inc. 384 Oyster Point Blvd., Suite 16 South San Francisco, CA 94080-1968 Tel: 1-800-553-4081 www.vitalwear.com

## SYSTEM SPECIFICATIONS:

PHYSICAL	
size (approximate)	8" diameter x 13.5" height (20.3 cm x 34.3 cm)
weight (dry)	5 lbs (2.27 kg)
fluid tether	8' (2.44 m)
CONTROL SYSTEM	
type	manual user control
THERMAL SYSTEM	
range	40°- 105° F (4.4°- 40.6° C)
contrast therapy transition time	1 min typical from extremes in temperature range
operating period	2-6 hours typical continuous use
thermal cutout	116.6° F (47° C)
CIRCULATING SYSTEM	
CIRCULATING SYSTEM	1.1 gal (4.2 L)
	1.1 gal (4.2 L) ice water
reservoir capacity	
reservoir capacity reservoir fluid	ice water
reservoir capacity reservoir fluid flow rate thru VW	ice water 5.5 gph (21 Lph)
reservoir capacity reservoir fluid flow rate thru VW maximum pressure	ice water 5.5 gph (21 Lph)
reservoir capacity reservoir fluid flow rate thru VW maximum pressure ELECTRICAL SYSTEM:	ice water 5.5 gph (21 Lph) 18 psi
reservoir capacity reservoir fluid flow rate thru VW maximum pressure ELECTRICAL SYSTEM: voltage	ice water 5.5 gph (21 Lph) 18 psi 120 V <sub>ac</sub> , 60 Hz
reservoir capacity reservoir fluid flow rate thru VW maximum pressure ELECTRICAL SYSTEM: voltage power	ice water 5.5 gph (21 Lph) 18 psi 120 V <sub>ac</sub> , 60 Hz 300 W
reservoir capacity reservoir fluid flow rate thru VW maximum pressure ELECTRICAL SYSTEM: voltage power leakage current	ice water 5.5 gph (21 Lph) 18 psi 120 V <sub>ac</sub> , 60 Hz 300 W
reservoir capacity reservoir fluid flow rate thru VW maximum pressure ELECTRICAL SYSTEM: voltage power leakage current CLASSIFICATION:	ice water 5.5 gph (21 Lph) 18 psi 120 V <sub>ac</sub> , 60 Hz 300 W under 300 μA

00630 Rev. E

For more information or to order VitalWrap Systems, wraps or accessories, contact our sales representatives at 1-800-553-4081 or visit www.vitalwear.com.